

Super-AMP Universal SINGLE TUBE 2.5X- RT/PCR PREMIX
(Mg final Conc. 3.5 mM) { Suitable for PCR/Real Time PCR}

Cat No.110206

DATA-SHEET

Ready to use pre-mixes are the innovation for convenience of your routine PCR as you just have to add your template RNA/DNA besides your primers & put the tubes for amplifications. These pre-mixes are optimized for most of the PCR templates. However amplifications may vary for different templates, hence Mg Concentration needs to be adjusted for further optimization. We are proud to Introduce our Super-AMP Universal Single Tube 2.5 X RT/PCR premix which will not only save you lot of time & money but will also help in controlling cross contaminations as cDNA synthesis & subsequent PCR can be done in the same tube without opening the tubes till the final amplification process.

Super-AMP Universal Single Tube 2.5X RT/PCR PRE-MIX is an optimized, stabilized, ready to use mixture of **blended M-Mulv**, Hot Start *Taq DNA polymerase*, RT/PCR buffer, MgCl₂ and dNTP's besides certain additives for enhanced amplifications. This 2.5X Single tube RT/PCR pre-mix performs cDNA synthesis & then subsequent PCR from a RNA template. The mixture is suitable for cDNA synthesis & subsequent amplification of most of the RNA templates for Conventional Single tube RT/PCR or Real Time Single tube RT/PCR. This premix is tested for absence of endo & exo nucleases and is **also tested for amplification of single gene copy**.

CONTENTS OF THE PACK :

Super-AMP Universal Single tube 2.5X RT/PCR-Premix 1.0 ml (for 50 reactions of 50µl)
Mol. Grade water 1.0 ml - 1 Vial

NOTE: It is not advisable to Freeze thaw the premix time & again. If the premix is to be used in several batches then it is advisable to prepare aliquots & store at -20 °C.

Mol grade water : Water is supplied in 1.0 ml quantity, which is RNAse and DNAse free. Advisable to aliquot the Mol. Grade water.

SPECIFICATIONS : Supplied as 2.5X in liquid form so as to ensure that you can add your primers & Template RNA volumes, to desired levels.

Storage Conditions : Recommended storage temperature is -20 °C.

DILUTION PROTOCOL for Super-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix:

Component	50 µl reaction Volume		25 µl reaction volume	
	Volume	Final Conc.	Volume	Final Conc.
2.5 X Super-Amp RT/PCR Premix	20 µl	1 X	10 µl	1 X
Forward primer	Variable	0.1-1 µM	Variable	0.1-1 µM
Reverse Primer	Variable	0.2-1 µM	Variable	0.2-1 µM
Template RNA	Variable	10pg-1 µg	Variable	10pg-1 µg
Sterile Deionized water	Up to 50 µl	----	Up to 25 µl	----

NOTE: -- Please use Gene specific primers for cDNA synthesis as well as for PCR.
-- It is advisable to use 20µl of extracted template RNA in case low yield is expected for 50µl reaction volumes. In case of High yields the volumes can be adjusted as per the yield of RNA.



Manufactured and Marketed By :

GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : RTP604/001

www.genome-diagnostics.com

www.genomediagnosics.in

Tel.: 00 91 01792 234285,, 234286

www.genomediagnosics.co.in

www.diagnosticsgenome.com

Suggested cycling protocol :

cDNA Synthesis	50°C	15 minutes	
Initial Denaturation	95°C	5 minutes	
Cycling	95°C	20 secs	40 – 45 cycles
	Annealing	30 secs	
	72°C extension	30 secs	
Final Extension	72°C	5 minutes	

- Note:**
- In case of long or highly structured templates the initial cDNA synthesis step can be extended to 30 minutes.
 - In case of gene specific primers if the Tm of the primers is above 55°C then the cDNA synthesis can be performed at 55°C also for 15 minutes.
 - In case of low copy number the Cycling can be extended to 45 Cycles.


Packing : The above Super-AMP Universal Single Tube 2.5X RT/PCR PRE-MIX is available in convenient small packs as well as large packs for bulk users.

The Super-AMP Universal Single Tube 2.5X RT/PCR PRE-MIX is supplied with an additional vial of DNase free water. Also RT/PCR Pre-mixes with various concentrations of MgCl₂ are available as mentioned below. Generally 2.5mM RT/PCR Pre-mix is considered as universal.

Pack size available :

Cat. No.	Product	Pack Size
110201	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 50 reactions of 50µl or 100 reactions of 25µl 1.5 mM final Mg at 1 X Conc.	1ml (50 rxns 50µl)
110202	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 100 reactions of 50µl or 200 reactions of 25µl 1.5 mM final Mg at 1 X Conc.	1ml x 2 (100 rxns 50µl)
110203	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 50 reactions of 50µl or 100 reactions of 25µl 2.5 mM final Mg at 1 X Conc.	1ml (50 rxns 50µl)
110204	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 100 reactions of 50µl or 200 reactions of 25µl 2.5 mM final Mg at 1 X Conc.	1ml x 2 (100 rxns 50µl)
110205	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 500 rxns of 50µl or 1000 reactions of 25µl 2.5 mM final Mg at 1 X Conc.	1ml x 10 (500 rxns 50µl)
110206	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 50 reactions of 50µl or 100 reactions of 25µl 3.5 mM final Mg at 1 X Conc.	1ml (50 rxns 50µl)
110207	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 100 reactions of 50µl or 200 reactions of 25µl 3.5 mM final Mg at 1 X Conc.	1ml x 2 (100 rxns 50µl)
110208	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 50 reactions of 50µl or 100 reactions of 25µl 4 mM final Mg at 1 X Conc.	1ml (50 rxns 50µl)
110209	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 100 reactions of 50µl or 200 reactions of 25µl 4 mM final Mg at 1 X Conc.	1ml x 2 (100 rxns 50µl)
110210	SUPER-AMP Universal SINGLE TUBE 2.5 X RT/PCR-Premix: 100 reactions of 50µl or 200 reactions of 25µl 5 mM final Mg at 1 X Conc.	1ml x 2 (100 rxns 50µl)

NOTE : For research use only. Not for use in Diagnostic or therapeutic procedures. Genome Diagnostics Pvt. Ltd., shall not in any event be liable for incidental or special damage of any kind resulting from any use except on applications.

	<p>Manufactured and Marketed By :</p> <p>GENOME DIAGNOSTICS PVT. LTD., (An ISO 13485:2003 & ISO 9001:2000 Certified Company) Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road, Near sector 4, Parwanoo, Dist. Solan. H.P. Email: genome24@rediffmail.com Version : RTP604/001</p> <p style="text-align: right;">Tel.: 00 91 01792 234285,, 234286</p> <p> www.genome-diagnostics.com www.genomediagnosics.co.in www.genomediagnosics.in www.diagnosticsgenome.com </p>
---	---

GENOME DIAGNOSTICS PVT. LTD.
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)



Manufactured and Marketed By :

GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : RTP604/001

Tel.: 00 91 01792 234285,, 234286

www.genome-diagnostics.com

www.genomediagnostics.co.in

www.genomediagnostics.in

www.diagnosticsgenome.com