

COMPLETE TEACHING KIT FOR PCR AMPLIFICATION
25 rxns of 50 µl OR 50 rxns of 25 µl

Cat No.110231

DATA-SHEET

Polymerase Chain Reaction (**PCR**) is the most powerful technique in molecular biology for amplification of DNA. By PCR one Copy of DNA can be amplified to billions of copies & can then be detected by various methods. We here provide a complete teaching kit for PCR amplification. This kit is specially designed for the Lecturers/Teachers to teach the students about the basic PCR amplification & beginners in the area of Molecular Biology. The kit provides an easy, optimized, PCR assay. The protocol described below is for 50µl reaction volume. In case 25µl reaction Volume is to be used then the Volumes can be just halved for the amplification reagents as shown in the table.

Contents of the Pack:

– Taq DNA Polymerase	65 U
– 10 X PCR Buffer	125 µl.
– Magnesium Chloride Sol.	0.5 ml
– dNTP Mix 40 mM	25 µl
– Mol. Grade Water	0.5 ml
– Lamda Primer Forward	50 µl
– Lamda Primer Reverse	50 µl
– Lamda DNA High Conc.	220 µl
– Lamda DNA Low Conc.	220 µl
– 100 bp DNA Mol Weight Marker	30 µl
– Gel Loading Dye	150 µl
– Mineral Oil	0.7 ml.
– Ethidium Bromide	50 µl
– 50 X TAE Buffer	50 ml.
– 0.2 ml PCR Tubes	25 Nos.
– Agarose 12.5 gms.	1 bottle.
– Pack Insert	1 Nos.

Storage of the Kit: Box No. 1 to be stored at -20°C. & Box No 3 to be stored at Room Temp alongwith Box No. 2.



Manufactured and Marketed By :

GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : TCH04/001

Tel.: 00 91 01792 234285,, 234286

www.genome-diagnostics.com

www.genomediagnostics.co.in

www.genomediagnostics.in

www.diagnosticsgenome.com

Preparation of Reagents:

All reagents are ready to use except for the following & their preparation is described below.

1 X TAE Buffer: To prepare 1 X TAE Buffer from 50 X dilute the provided 50 X TAE Buffer 1:50 with deionized water. i.e. 1ml of 50 X TAE buffer into 49 ml of Deionized Water.

Precautions:

- The UV in the laminar flow be put on before the usage of laminar flows at least for one hour.
- Turn off the UV at least 5 minutes before the usage of the Laminar flow.
- Put the laminar blower and Fluorescent lamp on.
- The work area should be set in a manner so that there is no hindrance in working.
- Thaw all components thoroughly at room temperature before starting.
- When thawed, mix the components and centrifuge briefly.
- Work quickly on ice or in the Cooling Block.
- Use pipette tips with filters only.
- Always use disposable powder-free gloves
- Avoid cross contamination
- Spilled droplet must be cleaned with Hypo
- Always discard all used tips etc. in the Hypo.

Protocol

Preparation for amplification:

	50 µl reaction Volume		25 µl reaction Volume	
	For 1 reaction	For 5 reactions	For 1 reaction	For 5 reactions
10 X PCR buffer	5.0 µl	25	2.5 µl	12.5
dNTP's Mix	1.0 µl	5 µl	0.5 µl	2.5 µl
MgCl ₂	4.0 µl	20 µl	2.0 µl	10 µl
Lambda Primer Forward	2.0 µl	10 µl	1.0 µl	5 µl
Lambda Primer Reverse	2.0 µl	10 µl	1.0 µl	5 µl
Taq DNA polymerase	0.5 µl	2.5 µl	0.25 µl	1.25 µl
Water	15.5 µl	77.5 µl	7.75 µl	38.75 µl
Total Volume	30 µl	150 µl	15 µl	75 µl



Manufactured and Marketed By :

GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : TCH04/001

Tel:: 00 91 01792 234285,, 234286

www.genome-diagnostics.com

www.genomediagnostics.co.in

www.genomediagnostics.in


www.diagnosticsgenome.com

Protocol for 50 µl reaction Volume


-- Place the 0.2 ml PCR tubes provided in this kit in the PCR rack (Not provided in kit) & dispense 30 µl of the above mixture in each PCR tube.

-- Dispense 20 µl of Low Conc. Lamda DNA (in two tubes) & High Conc. Lamda DNA (in two tubes) & 20 µl of Mol Grade water in one tube as described below.

20 µl of Low Conc. Lamda DNA



20 µl of High Conc. Lamda DNA



20 µl of Mol. Grade Water.



NOTE: In case of Hot Lid thermal Cyclers there is no need for adding Mineral oil. In case the Thermal Cycler is not having Hot Lid facility then please overlay 25 µl of Mineral oil in each PCR tube by adding the Mineral oil through the sides of the PCR tube as depicted below.



Close the cap of the PCR tubes & load these 0.2ml tubes in the PCR machine and set the cycling conditions as follows:

Initial denaturation : 95° C for 5 min

Cycling:

Denaturation : 95° C for 30 sec
Annealing : 55° C for 30 sec
Extension : 72° C for 30 sec } For 40 cycles

Final extension : 72° C for 5 min



Manufactured and Marketed By :

GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : TCH04/001

Tel.: 00 91 01792 234285,, 234286

www.genome-diagnostics.com


www.genomediagnostics.co.in

www.genomediagnostics.in

www.diagnosticsgenome.com

After PCR, perform Agarose gel electrophoresis as described below:

1. Add 1 ml of 50 X TAE buffer to 49 ml of distilled water in order to obtain 1X TAE buffer as described in Preparation of reagents.
2. Add 600 mg of Agarose (provided one pouch) to 50 ml of freshly prepared 1X TAE buffer in a conical Glass flask.
3. Boil the mixture to obtain a clear Agarose solution. Alternatively it can be boiled in a microwave oven as well. As soon as the Sol is clear stop boiling.
4. Add 2 μ l of ethidium bromide solution to the agarose solution & mix gently.
5. Pour this solution into a clean gel-casting tray, which has already been sealed from two open sides for the casting of Agarose gel. Place the combs for the formation of the wells in appropriate position immediately after pouring the agarose solution.
6. Allow the mixture to cool for 20 – 30 min to form a gel.
7. After the formation of gel, remove the comb and the edge sealing & transfer the gel to the Electrophoretic gel tank.
8. Fill the electrophoresis gel tank with 1X TAE buffer prepared as mentioned in Preparation of reagents. (Total volume of 1 X TAE to be poured depends on the capacity of the Gel Tank, which could be between 750 ml to 1.5 liters. Please refer to the manual of your Gel Electrophoresis system for confirming the tank volumes)
9. Add 10 μ l of the gel loading dye (provided) to the PCR mixture in each tube, which was obtained after PCR amplification. Mix them thoroughly. (In case Mineral oil has been overlaid then remove the mineral oil first with the help of a micropipette & then add & mix the Gel loading Dye in the PCR tubes.
10. Load 6 μ l of the marker (DNA Ladder – provided) in the first well of the Agarose gel.
11. Load 20 μ l of the PCR product mixed with gel loading dye from above tubes from the second well in the Agarose gel onwards as depicted below.



12. Apply 70 – 80V for 45 minutes-1 hour.
13. After the electrophoresis, transfer the gel on to the Transilluminator.
14. View the gel at 260/280 nm.
15. You must be able to observe the bands ladder in the first lane, and a faint band of 277 bp in lane no. 2 & 3, & a bright band of 277 bp in lane no 4 & 5 & no band in lane no 6 which confirms the success of PCR procedure as shown below.



Manufactured and Marketed By :

GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : TCH04/001


Tel:: 00 91 01792 234285,, 234286

www.genome-diagnostics.com

www.genomediagnostics.co.in

www.genomediagnostics.in

www.diagnosticsgenome.com



Special Note:

1. The TAE Buffer once prepared can be Used for at least 5 subsequent runs. Please store the same in a bottle after the run.
2. The Lamda DNA & primers have been provided for the Convenience of Usage. However if you wish to your use your own Primers and DNA then it can simply be replaced by the Lamda Primers and DNA and Lamda controls.

NOTE : For research use only. Not for use in Diagnostic or therapeutic procedures. Genome Diagnostics Pvt. Ltd., shall not in any event be liable for incidental or special damage of any kind resulting from any use except on applications.



Manufactured and Marketed By :

GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : TCH04/001

Tel:: 00 91 01792 234285,, 234286

www.genome-diagnostics.com

www.genomediagnosics.co.in

www.genomediagnosics.in

www.diagnosticsgenome.com