

## HOT START *Taq* DNA Polymerase ( 250 Units)

Cat No. 110042

### DATA-SHEET

Hot Start *Taq* DNA polymerase is an optimized mixture of *Taq* DNA polymerase and anti-*Taq* DNA polymerase monoclonal antibodies. Antibodies block polymerase activity during set-up of the PCR reactions at ambient temperature (20-22°C). The inhibition of *Taq* DNA Polymerase is reversed when the temperature is above 70°C. The PCR products obtained with Hot Start *Taq* are free from non-specific products and primer-dimers.

### UNIT DEFINITION :

One unit of activity is defined as the amount of enzyme required to incorporate 10 nmoles of dNTP's into acid-insoluble DNA fraction in 30 minutes at 72°C.

**CONCENTRATION SUPPLIED :** Enzyme is supplied at a Concentration of 5U/μl along with one vial of 10X buffer, one vial of additional MgCl<sub>2</sub> & one vial of DNA's free water.

**STORAGE BUFFER :** 10mM Tris-HCL (pH 7.0), 0.1mM EDTA, 50mM KCl, 50% glycerol.

### CONTENTS OF THE PACK :

250 units of Hot Start <i>Taq</i> DNA polymerase	50 μl	- 1 Vial
10 X Reaction buffer	0.5 ml	- 1 Vial
25 mM MgCl <sub>2</sub>	1.5 ml	- 1 Vial
Mol. Grade water	1.5 ml	- 1 Vials

### 10X REACTION BUFFER COMPLETE INCLUDING MAGNESIUM (2.5 mM final 1x Conc.)

670mM Tris-HCL (pH 8.8), 160mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 0.1% Tween-20, 25mM MgCl<sub>2</sub>.

**MgCl<sub>2</sub> Solution :** 1.5 ml of 25mM MgCl<sub>2</sub> Solution (In case higher Conc. of Mg is to be used)

**Mol grade water:** Water is supplied in 1.5 ml X 1 vials, which is RNase and DNase free. It is advisable to aliquot the Mol. Grade water.

**Storage Conditions :** Recommended storage temperature is – 20 °C.

### AVIALABLE PACK SIZE :

Cat. No.	PRODUCT	Pack Size
110041	Hot Start <i>Taq</i> DNA Polymerase	100 U
110042	Hot Start <i>Taq</i> DNA Polymerase	250 U
110043	Hot Start <i>Taq</i> DNA Polymerase	500 U (250 U x 2)
110044	Hot Start <i>Taq</i> DNA Polymerase	1000 U (250 U x 4)
110045	Hot Start <i>Taq</i> DNA Polymerase	5000 U (1000 U x 5)

**NOTE :** For research use only. Not for use in Diagnostic or therapeutic procedures.

Genome Diagnostics Pvt. Ltd. shall not in any event be liable for incidental or special damage of any kind resulting from any use except on applications.



Manufactured and Marketed By :

**GENOME DIAGNOSTICS PVT. LTD.,**  
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,  
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : M04/001

[www.genome-diagnostics.com](http://www.genome-diagnostics.com)

[www.genomediagnosics.in](http://www.genomediagnosics.in)

Tel.: 00 91 01792 234285,, 234286

[www.genomediagnosics.co.in](http://www.genomediagnosics.co.in)

[www.diagnosticsgenome.com](http://www.diagnosticsgenome.com)