

110253	dNTP mix 100 mM: (dATP, dCTP, dGTP and dTTP each, at a concentration of 25mM)	1ml x 5
110254	dNTP mix 100 mM: (dATP, dCTP, dGTP and dTTP each, at a concentration of 25mM)	1ml x 10

INDIVIDUAL dNTP's solutions(40 mM & 100 mM)

Cat. No.	PRODUCT	Pack Size
110261	dATP 40mM	500 µ l
110262	dATP 40mM	1ml
110263	dATP 100mM	500 µ l
110264	dATP 100mM	1ml
110265	dCTP 40mM	500 µ l
110266	dCTP 40mM	1ml
110267	dCTP 100mM	500 µ l
110268	dCTP 100mM	1ml
110269	dGTP 40mM	500 µ l
110270	dGTP 40mM	1ml
110271	dGTP 100mM	500 µ l
110272	dGTP 100mM	1ml
110273	dTTP 40mM	500 µ l
110274	dTTP 40mM	1ml
110275	dTTP 100mM	500 µ l
110276	dTTP 100mM	1ml

Set of 4 dNTP's (dATP, dGTP, dCTP, dTTP)

Cat. No.	PRODUCT	Pack Size
110281	dATP, dCTP, dGTP, dTTP 40mM each	200 µ l x 4
110282	dATP, dCTP, dGTP, dTTP 40mM each	1ml x 4
110283	dATP, dCTP, dGTP, dTTP 100mM each	200 µ l x 4
110284	dATP, dCTP, dGTP, dTTP 100mM each	1ml x 4

dUTP

dUTP may be used instead of dTTP in the PCR and RT-PCR protocols to prevent carryover from previous amplifications. The substitution of dUTP for dTTP in PCR results in uracil-containing PCR products that are suitable for most standard applications. The Enzyme Uracil-DNA glycosylase, UDG can be added to a PCR premix to excise uracil from any contaminating PCR product, thereby preventing false positives.

Cat. No.	PRODUCT	Pack Size
110291	dUTP 40mM	500 µ l
110292	dUTP 40mM	1ml
110293	dUTP 100mM	500 µ l
110294	dUTP 100mM	1ml

POWDERED dNTP's & dUTP

Cat. No.	PRODUCT	Pack Size
110301	dATP powder	20 mg
110302	dATP powder	100 mg
110303	dATP powder	1 gm
110304	dCTP powder	20 mg

Manufactured and Marketed By :



GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : P104/001

www.genome-diagnostics.com

www.genomediagnostics.in

Tel: 00 91 01792 234285,, 234286

www.genomediagnostics.co.in

www.diagnosticsgenome.com

GENOME DIAGNOSTICS PVT. LTD.
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

110305	dCTP powder	100 mg
110306	dCTP powder	1 gm
110307	dGTP powder	20 mg
110308	dGTP powder	100 mg
110309	dGTP powder	1 gm
110310	dTTP powder	20 mg
110311	dTTP powder	100 mg
110312	dTTP powder	1 gm
110313	dUTP powder	20 mg
110314	dUTP powder	100 mg
110315	dUTP powder	1 gm

NOTE : For research use only. Not for use in Diagnostic or therapeutic procedures.
Genome Diagnostics Pvt. Ltd., shall not in any event be liable for incidental or special damage of any kind resulting from any use except on applications.

Manufactured and Marketed By :



GENOME DIAGNOSTICS PVT. LTD.,
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: genome24@rediffmail.com

Version : P104/001

www.genome-diagnostics.com

www.genomediagnostics.in

Tel: 00 91 01792 234285,, 234286

www.genomediagnostics.co.in

www.diagnosticsgenome.com