

**Proteinase K 20mg/ml**  
5ml

Cat No.110322

**DATA-SHEET**

Proteinase K is a Broad-Spectrum serine protease for general digestion of proteins in Biological samples. Enzyme is free from RNase and DNase activities. Proteinase K is stable in urea and SDS. Recommended working concentration of Proteinase K is 10-20 µg/ml for majority of applications.

**CONTENTS OF THE PACK :**

20mg/ml of Proteinase K                      1ml - 5 Vials

**STORAGE :** Recommended storage is at -20°C.

**AVAILABILITY :** Both in Liquid ready to use solution and in powder form.

Cat. No.	PRODUCT	Pack Size
110321	Proteinase K liquid 20mg/ml	1ml
110322	Proteinase K liquid 20mg/ml	1 ml x 5
110323	Proteinase K liquid 20mg/ml	1 ml x 10
110331	Proteinase K Powder	20 mg
110332	Proteinase K Powder	100 mg
110333	Proteinase K Powder	1 gm

**NOTE :** For research use only. Not for use in Diagnostic or therapeutic procedures.  
Genome Diagnostics Pvt. Ltd., shall not in any event be liable for incidental or special damage of any kind resulting from any use except on applications.



Manufactured and Marketed By :

**GENOME DIAGNOSTICS PVT. LTD.,**  
(An ISO 13485:2003 & ISO 9001:2000 Certified Company)

Up Mohal Naryal, Khasra No. 427, Opp. Divya Packers, Old Timber Depot Road,  
Near sector 4, Parwanoo, Dist. Solan. H.P.

Email: [genome24@rediffmail.com](mailto:genome24@rediffmail.com)

Tel: 00 91 01792 234285,, 234286

[www.genome-diagnostics.com](http://www.genome-diagnostics.com)

[www.genomediagnostics.co.in](http://www.genomediagnostics.co.in)

[www.genomediagnostics.in](http://www.genomediagnostics.in)

[www.diagnosticsgenome.com](http://www.diagnosticsgenome.com)